

Product Catalog

July 2019

Nucleic acid analogs for: DNA sequencing Diagnostics SNP detection Synthetic Biology

www.firebirdbio.com

Contents

	Reversible terminators DNA sequencing and SNP identification	4
s, LLC	AEGIS TM Artificially Expanded Genetic Information Systems Orthogonality, cleanliness, and uniformity in nucleic acid analysis	10
	SAMRS TM Self-Avoiding Molecular Recognition Systems Reducing primer dimers in nucleic acid amplification	19
N SC	Biversals Novel bases for binding to degenerate DNA sequences	22
CULA	SNAP-2 Identifying DNA and RNA molecules when the exact sequence is not known	24
OLI	Ordering	25
	References	26

Firebird reagents and the origins of life

This catalog presents reagents and technologies to aid in the development of new biomedical innovations. Several of the products presented here were developed with the goal of better understanding the features of DNA and RNA that make them the carriers of genetic information for life on earth (Karalkar and Benner, 2018).

The figure shows the crystal structure of a double helix built from an expanded genetic alphabet of eight hachimoji building blocks, G (green), A (red), C (dark blue), T (yellow), B (aqua) P (blue), S (violet), and Z (orange). Adopted from Hoshika et al., 2019.

Reagents in this catalog are sold for research use only. Please inquire for IVD use.

Reversible terminators

Firebird has introduced the 3'-ONH₂ reversible terminator as an alternative to the larger 3'-OCH₂N₃ group (Hutter et al., 2010). The 3'-ONH₂ group is accepted by a variety of polymerases. After incorporation, further primer extension is blocked. The 3'-ONH₂ group is cleaved with buffered aqueous sodium nitrite to regenerate a standard 3'-OH.



3'-ONH₂ reversible terminators can be used in DNA sequencing (Hutter et al., 2010), Oligonucleotide synthesis (Jensen and Davis, 2018, Sarac and Hollenstein, 2019) and SNP analysis (Chen et al., 2010).

Polymerase variants have been developed that incorporate these terminators with improved efficiency over standard Taq DNA polymerase (Chen et al., 2010):

POL475-400	400 units	\$110.00
POL475-1000	1000 units	\$250.00

Reversible terminators: Which one?

As shown on the following pages, Firebird offers reversible terminators in several different formats:

The ONH2-170 series: Ready to use, untagged The ONH2-100 series: Oxime blocked, untagged The ONH2-DT series: Diol-linked tags The ONH2-DS series: Disulfide-linked tags

The 170 series ("ready to use") are crude products from a solid phase synthesis procedure that results in approximately 90% pure terminator triphosphates, which should be sufficient for most purposes. This series carries a free 3'-ONH₂ group and can be incorporated by various enzymes without further modification. The shelf life of this series is about one year (if stored properly, at or below -20°C).

The 100 series ("oxime blocked") are the acetoxime form of the 170 series and has a longer shelf life (many years, if stored properly). This series is HPLC purified to remove most of the impurities, albeit at an increased cost. But these oximes will not be incorporated by enzymes. For enzymatic incorporation the free 3'-ONH₂ must be regenerated before use, which can be achieved in situ with buffered aqueous methoxylamine.

The DT and DS series are reversible terminators carrying a cleavable linker attached to the nucleobase with a free amino group, allowing for the attachment of a tag (e.g. fluorescent dye or biotin). As in the 100 series, the 3'-ONH₂ group is protected as the acetoxime, for attachment of the tag and subsequent purification as well as increased shelf life. And as with the 100 series, this oxime needs to be deprotected with buffered aqueous methoxylamine before use. Upon incorporation of the terminator into an oligonucleotide the linker can be cleaved under mild conditions (dilute aqueous periodate for the DT series, phosphine or thiol reagents for the DS series) to remove the tag.

For a fee, Firebird scientists can attach a customer specified tag onto the primary amine; please inquire for pricing.

Reversible terminators: Ready to use, untagged

Reversible terminators with a free 3'-ONH₂ are ready to use and do not require deprotection or further treatment before use.









Cytosine	
CONH2-172	
5 µmoles	\$750.00
50 µmoles	\$5250.00



Adenine	
AONH2-173	
5 µmoles	\$750.00
50 μmoles	\$5250.00



Guanine	
GONH2-174	
5 µmoles	\$750.00
50 µmoles	\$5250.00

Please inquire for bulk pricing.

Reversible terminators: Oxime blocked, untagged

Reversible terminators are also available with the 3'-ONH₂ protected as the acetoxime, which are more suitable for long term storage than our ready-to-use reversible terminators. For enzymatic incorporation, the oxime must be deprotected to the free 3'-ONH₂ before use, which can be achieved in situ with buffered aqueous methoxylamine.



Please inquire for bulk pricing and availability.

Reversible terminators: diol-linked tags

Firebird offers reversible terminators with a 3'-ONH₂ (protected as the acetoxime) and a diol linker carrying a free amino group, to which a tag (fluorescent dye or other moiety) can be attached. The diol can be rapidly cleaved with dilute aqueous periodate.



For enzymatic incorporation, the oxime must be deprotected to the free 3'-ONH₂ before use, which can be achieved in situ with buffered aqueous methoxylamine.

Please inquire for bulk pricing.

Not available for sequencing applications on certain machines and in certain jurisdictions. Please inquire.

Reversible terminators: disulfide-linked tags

Firebird offers reversible terminators with a 3'-ONH₂ (protected as the acetoxime) and a disulfide linker carrying a free amino group, to which a tag (fluorescent dye or other moiety) can be attached. The disulfide can be rapidly cleaved with phosphine or thiol reagents.





For enzymatic incorporation, the oxime must be deprotected to the free 3'-ONH₂ before use, which can be achieved in situ with buffered aqueous methoxylamine.

Please inquire for bulk pricing.

AEGISTM

Artificially Expanded Genetic Information Systems

Firebird has created different heterocycles in order to implement additional hydrogen bonding patterns for the Z:P and S:B base pairs as shown below. These base pairs are "orthogonal" to C:G and A:T (Benner, 2004; Hoshika et al., 2019; Sefah et al., 2014; Yang et al., 2006).





Using AEGIS: For ultra-clean nested PCR

The Z:P pair is retained during PCR amplification, even when present as consecutive base pairs (Yang et al., 2011).



AEGIS-nested primers suppress noise in multiplex PCR (Yang et al., 2010) and support 22-plex PCR to detect RNA viruses (Glushakova et al., 2015a).

"Analyte specific" chimeric primers in low concentration



Using AEGIS: Assembling long DNA constructs

Adding AEGIS nucleotides to single-stranded oligos increases their information density. This allows clean and rapid hybridization of single stranded DNA unobstructed by hairpins, wandering strands, and non-canonical interactions. Then, using Firebird's "transliteration" technology, the AEGIS nucleotides are cleanly replaced by standard nucleotides giving an entirely natural gene (Bradley and Benner, 2014; Merritt et al., 2014).

An example of this assembly is shown below: S and B residues in the tails of the primers guide specific alignment of the unique sequences at the ends of the primers. The SB basepair is converted to a TA basepair during PCR amplification.



Using AEGIS: Molecular beacons

Adding AEGIS nucleotides to the stems of molecular beacons prevents their being opened by adventitious DNA and RNA in a complex biological sample. This allows beacons to deliver signals in complex biological samples with significantly lower noise (Sheng et al., 2008).





Firebird offers custom-synthesized beacons with a wide range of fluorescent dyes and quenchers. Please inquire.

AEGIS: Increasing Aptamer diversity

The dP and dZ bases increase the chemical repertoire of nuleic acids for aptamer selections (Biondi et al., 2016; Sefah et al., 2014; Zhang et al., 2015, 2016). Of special note is the NO_2 group present on dZ.

The figure below shows one application of AEGIS for the selection of aptamers specific to a cell surface protein on a cancer cell (Zhang et al., 2016).



AEGIS: dZ and dP Phosphoramidites

dZ Phosphoramidite



Hydrogen bonding pattern: Z Sugar: 2'-Deoxyribose Heterocycle: Nitropyridine Linear Formula C₅₀H₅₇N₆O₁₂P Mol Weight 964.99

Cat. No. dZ-PA-101 100 mg \$ 720.00 1 gram \$5760.00

dP Phosphoramidite



Hydrogen bonding pattern: P Sugar: 2'-Deoxyribose Heterocycle: Imidazotriazine Linear Formula C₄₉H₆₅N₈O₇P Mol Weight 909.08

Cat. No. dP-PA-102 100 mg \$ 850.00 1 gram \$6800.00

(Yang et al., 2006) Please inquire for bulk pricing

Oligonucleotides are available that incorporate dZ and dP. Ribonucleoside derivatives are also available.

AEGIS: dZ and dP Triphosphates

dZ Triphosphate



Hydrogen bonding pattern: Z Sugar: 2'-deoxyribose Heterocycle: Nitropyridine Linear Formula C₁₀H₁₆N₃O₁₅P₃ Mol Weight 511.17

Cat. No. dZTP-101 1 µmole \$ 400.00 5 µmoles \$ 1600.00

dP Triphosphate



Hydrogen bonding pattern: P Sugar: 2'-deoxyribose Heterocycle: Imidazotriazine Linear Formula C₁₀H₁₆N₅O₁₃P₃ Mol Weight 493.15

Cat. No. dPT	P-20	01
1 μmole	\$	420.00
5 μmoles	\$	1680.00

(Yang et al., 2007) Please inquire for bulk pricing.

Ribonucleoside triphosphates are also available. Please inquire.

AEGIS: dS and dB Phosphoramidites

dS Phosphoramidite



Hydrogen bonding pattern: S Sugar: 2'-Deoxyribose Heterocycle: Pseudo-C Linear Formula C₄₇H₅₄N₅O₈P Mol Weight 847.93

Cat. No. dS-PA-104S 100 mg \$ 720.00 1 gram \$5760.00

dB Phosphoramidite



Hydrogen bonding pattern: B Sugar: 2'-Deoxyribose Heterocycle: Purine Linear Formula C₄₉H₆₈N₇O₇P Mol Weight 898.08

Cat. No. dB-PA-103P 100 mg \$ 600.00 1 gram \$4800.00

Please inquire for bulk pricing.

Oligonucleotides are available that incorporate dS and dB. Ribonucleoside derivatives are also available.

AEGIS: dS and dB Triphosphates

dS Triphosphate



Hydrogen bonding pattern: S Sugar: 2'-deoxyribose Heterocycle: Pseudo-C Linear Formula $C_{10}H_{18}N_3O_{13}P_3$ Mol Weight 481.18

Cat. No. dSTP-401S 1 μmole \$ 400.00 5 μmoles \$ 1600.00

dB Triphosphate



Hydrogen bonding pattern: B Sugar: 2'-deoxyribose Heterocycle: Purine Linear Formula C₁₀H₁₆N₅O₁₃P₃ Mol Weight 507.18

Cat. No. dB7	P-301P	
l µmole	\$ 400.00	
5 µmoles	\$ 1600.00	

Please inquire for bulk pricing.

Ribonucleoside triphosphates are also available. Please inquire.

SAMRSTM

Self Avoiding Molecular Recognition System

Nucleic acid amplification reactions suffer from the formation of off-target amplification products. This problem increases with the number of primers in the reaction. It mainly occurs because of primer-primer interactions. Thus, selectively removing hydrogen bonding units from the primer bases results in *self-avoiding* DNA primers that can be added to enable and improve multiplexed PCR (Hoshika et al., 2010).

As shown below, SAMRS bases, indicated with a *, can base pair with standard bases (in the target or amplicon) but not with other SAMRS bases.



SAMRS: design rules

General rules for placement of SAMRS bases into primers used in PCR and other amplification techniques. Upper case letters indicate the standard base, lower case letters indicate the SAMRS base.

- 1. Recommended lengths are 20-35 nts. SAMRS bases should be utilized in the first 4, (up to 8), positions at the 3' end of the oligonucleotide <u>but not in the very first 3' base</u>.
- 2. Between two and four SAMRS bases should be used per oligonucleotide, with two or three SAMRS modifications being preferred.
- 3. The SAMRS t base has not demonstrated as much reduction in primer dimer formation so it is preferable to substitute a, g or c rather than t when given a choice.
- 4. Avoid using SAMRS bases in a string of three and four consecutive identical bases. For example: ggg, ccc, aaa, ttt, gggg, cccc, aaaa, and tttt are to be avoided.
- 5. Separate SAMRS bases with standard bases, for example: gGg, cCc, aAa, and tTt.
- 6. If the 3'-end of an oligo is standard A or T, the second base should also be a standard base, followed by two or three SAMRS modifications.

(Benner et al., 2015; Glushakova et al., 2015a, 2015b; Hoshika et al., 2010; Sharma et al., 2014; Yang et al., 2015 and unpublished data)

As with all primer designs, Firebird makes no guarantee as to the performance of oligonucleotides made by following these rules.

SAMRS: Phosphoramidites

t* phosphoramidite 1 g \$ 1050 10 g \$ 8220

a* phosphoramidite 1 g \$ 655 10 g \$ 5340

c* phosphoramidite 1 g \$ 394 10 g \$3150

g* phosphoramidite 1 g \$ 110 10 g \$ 880

Oligonucleotides are available that incorporate the SAMRS bases.

SAMRS can be combined with AEGIS in oligos to get highly multiplexed and highly clean PCR, as well as isothermal amplification. Please inquire.

BiversalTM nucleotides

The divergence present in pathogen genomes makes it difficult or impossible to design a single probe or primer that is the exact complement of all possible target sequences. In practice this is handled using degenerate or "universal bases" (e.g. inosine).

Firebird scientists have created two "biversals" TM. Each exists in two tautomeric forms that provide alternate hydrogen bonding patterns. This allows the pyrimidine biversal to pair with either guanine or adenine and the purine biversal to pair with either thymine or cytosine. These bases provide a new alternative to the existing methods of handling divergent sequences for both PCR amplification and SNP detection (Yang et al., 2018).

Tautomeric forms of the Pyrimidine biversal:



Tautomeric forms of the Purine biversal:



Biversal: Phosphoramidites





Pyrimidine biversal phosphoramidite Linear Formula C₄₁H₅₀N₅O₈P Mol Weight 771.84

Cat No. Pyr-BiVer

100 mg	\$ 500
1 gram	\$ 4000

Purine biversal phosphoramidite Linear Formula C₅₀H₆₂N₇O₉P Mol Weight 936.04

Cat No. Pur-BiVer

100 mg	\$ 2850
1 gram	Inquire

Firebird sells oligonucleotides that contain the purine and pyrimidine biversals. Please inquire.

SNAP2[™] Oligonucleotides: Priming with high specificity

Specific DNA priming requires oligonucleotides with complementarity of ≥ 16 nucleotides but the ability to discriminate against single mismatches. Dynamic assembly of a primer on a template can provide this (Leal et al., 2006).





Controlled pore glass with aldehyde precursor SNAP2T_HJK001 100 mg \$600.00

Firebird sells phosphoramidites, oligos and libraries with SNAP2 ends. Please inquire.

SNAP2 primers are not to be confused with Snap-Tag® and other registered trademarks of New England BioLabs.

Purchasing, shipping and handling

Contact

Orders By email: orders@firebirdbio.com By telephone: (386) 418-0347 9:00 - 5:00 Eastern Time

<u>Technical inquiries</u> support@firebirdbio.com

Reagents and enzymes are shipped at the customer's expense, on dry ice as appropriate.

References

Benner, S.A. (2004). Chemistry. Redesigning genetics. Science 306, 625-626.

Benner, S.A., Kim, H.-J., Merritt, K.B., Yang, Z., McLendon, D.C., Hoshika, S., and Hutter, D. (2015). Next-generation DNA in pathogen detection, surveillance, and CLIA-waivable diagnostics. SPIE Digital Library 9490, 94900K-94900K – 6.

Biondi, E., Lane, J.D., Das, D., Dasgupta, S., Piccirilli, J.A., Hoshika, S., Bradley, K.M., Krantz, B.A., and Benner, S.A. (2016). Laboratory evolution of artificially expanded DNA gives redesignable aptamers that target the toxic form of anthrax protective antigen. Nucleic Acids Res. 44, 9565–9577, PMCID: PMC5175368.

Bradley, K.M., and Benner, S.A. (2014). OligArch: A software tool to allow artificially expanded genetic information systems (AEGIS) to guide the autonomous self-assembly of long DNA constructs from multiple DNA single strands. Beilstein J Org Chem 10, 1826–1833, PMCID: PMC4142867.

Chen, F., Gaucher, E.A., Leal, N.A., Hutter, D., Havemann, S.A., Govindarajan, S., Ortlund, E.A., and Benner, S.A. (2010). Reconstructed evolutionary adaptive paths give polymerases accepting reversible terminators for sequencing and SNP detection. Proc. Natl. Acad. Sci. U.S.A. 107, 1948– 1953, PMCID: PMC2804741.

Glushakova, L.G., Bradley, A., Bradley, K.M., Alto, B.W., Hoshika, S., Hutter, D., Sharma, N., Yang, Z., Kim, M.-J., and Benner, S.A. (2015a). Highthroughput multiplexed xMAP Luminex array panel for detection of twenty two medically important mosquito-borne arboviruses based on innovations in synthetic biology. J. Virol. Methods 214, 60–74, PMCID: PMC4485418.

Glushakova, L.G., Sharma, N., Hoshika, S., Bradley, A.C., Bradley, K.M., Yang, Z., and Benner, S.A. (2015b). Detecting respiratory viral RNA using expanded genetic alphabets and self-avoiding DNA. Anal. Biochem. 489, 62– 72, PMCID: PMC4733849.

Hoshika, S., Chen, F., Leal, N.A., and Benner, S.A. (2010). Artificial genetic systems: self-avoiding DNA in PCR and multiplexed PCR. Angew. Chem. Int. Ed. Engl. 49, 5554–5557, PMCID: PMC6027612.

Hoshika, S., Leal, N.A., Kim, M.-J., Kim, M.-S., Karalkar, N.B., Kim, H.-J., Bates, A.M., Watkins, N.E., SantaLucia, H.A., Meyer, A.J., et al. (2019). Hachimoji DNA and RNA: A genetic system with eight building blocks. Science 363, 884–887, PMCID: PMC6413494.

References

- Hutter, D., Kim, M.-J., Karalkar, N., Leal, N.A., Chen, F., Guggenheim, E., Visalakshi, V., Olejnik, J., Gordon, S., and Benner, S.A. (2010). Labeled nucleoside triphosphates with reversibly terminating aminoalkoxyl groups. Nucleosides Nucleotides Nucleic Acids 29, 879–895, PMCID: PMC3858015.
- Jensen, M.A., and Davis, R.W. (2018). Template-Independent Enzymatic Oligonucleotide Synthesis (TiEOS): Its History, Prospects, and Challenges. Biochemistry 57, 1821–1832.
- Karalkar, N.B., and Benner, S.A. (2018). The challenge of synthetic biology. Synthetic Darwinism and the aperiodic crystal structure. Curr Opin Chem Biol 46, 188–195.
- Leal, N.A., Sukeda, M., and Benner, S.A. (2006). Dynamic assembly of primers on nucleic acid templates. Nucleic Acids Res. 34, 4702–4710, PMCID: PMC1635275.
- Merritt, K.K., Bradley, K.M., Hutter, D., Matsuura, M.F., Rowold, D.J., and Benner, S.A. (2014). Autonomous assembly of synthetic oligonucleotides built from an expanded DNA alphabet. Total synthesis of a gene encoding kanamycin resistance. Beilstein J Org Chem 10, 2348–2360, PMCID: PMC4222377.
- Sarac, I., and Hollenstein, M. (2019). Terminal Deoxynucleotidyl Transferase in the Synthesis and Modification of Nucleic Acids. Chembiochem 20, 860–871.
- Sefah, K., Yang, Z., Bradley, K.M., Hoshika, S., Jiménez, E., Zhang, L., Zhu, G., Shanker, S., Yu, F., Turek, D., et al. (2014). In vitro selection with artificial expanded genetic information systems. Proc. Natl. Acad. Sci. U.S.A. 111, 1449–1454, PMCID: PMC3910645.
- Sharma, N., Hoshika, S., Hutter, D., Bradley, K.M., and Benner, S.A. (2014). Recombinase-based isothermal amplification of nucleic acids with self-avoiding molecular recognition systems (SAMRS). Chembiochem 15, 2268–2274.
- Sheng, P., Yang, Z., Kim, Y., Wu, Y., Tan, W., and Benner, S.A. (2008). Design of a novel molecular beacon: modification of the stem with artificially genetic alphabet. Chem. Commun. (Camb.) 5128–5130, PMCID: PMC2763601.

References

- Yang, Z., Hutter, D., Sheng, P., Sismour, A.M., and Benner, S.A. (2006). Artificially expanded genetic information system: a new base pair with an alternative hydrogen bonding pattern. Nucleic Acids Res. 34, 6095– 6101, PMCID: PMC1635279.
- Yang, Z., Sismour, A.M., Sheng, P., Puskar, N.L., and Benner, S.A. (2007). Enzymatic incorporation of a third nucleobase pair. Nucleic Acids Res. 35, 4238–4249, PMCID: PMC1934989.
- Yang, Z., Chen, F., Chamberlin, S.G., and Benner, S.A. (2010). Expanded genetic alphabets in the polymerase chain reaction. Angew. Chem. Int. Ed. Engl. 49, 177–180, PMCID: PMC3155763.
- Yang, Z., Chen, F., Alvarado, J.B., and Benner, S.A. (2011). Amplification, mutation, and sequencing of a six-letter synthetic genetic system. J. Am. Chem. Soc. 133, 15105–15112, PMCID: PMC3427765.
- Yang, Z., McLendon, C., Hutter, D., Bradley, K.M., Hoshika, S., Frye, C.B., and Benner, S.A. (2015). Helicase-Dependent Isothermal Amplification of DNA and RNA by Using Self-Avoiding Molecular Recognition Systems. Chembiochem 16, 1365–1370, PMCID: PMC4489552.
- Yang, Z., Kim, H.-J., Le, J.T., McLendon, C., Bradley, K.M., Kim, M.-S., Hutter, D., Hoshika, S., Yaren, O., and Benner, S.A. (2018). Nucleoside analogs to manage sequence divergence in nucleic acid amplification and SNP detection. Nucleic Acids Res. 46, 5902–5910, PMCID: PMC6159519.
- Zhang, L., Yang, Z., Sefah, K., Bradley, K.M., Hoshika, S., Kim, M.-J., Kim, H.-J., Zhu, G., Jiménez, E., Cansiz, S., et al. (2015). Evolution of functional six-nucleotide DNA. J. Am. Chem. Soc. 137, 6734–6737, PMCID: PMC4500535.
- Zhang, L., Yang, Z., Le Trinh, T., Teng, I.-T., Wang, S., Bradley, K.M., Hoshika, S., Wu, Q., Cansiz, S., Rowold, D.J., et al. (2016). Aptamers against Cells Overexpressing Glypican 3 from Expanded Genetic Systems Combined with Cell Engineering and Laboratory Evolution. Angew. Chem. Int. Ed. Engl. 55, 12372–12375, PMCID: PMC5554412.